


**DR. 365 PHARMACEUTICAL & EQUIPMENTS**



Web: [www.dr365pharmaceutival.com](http://www.dr365pharmaceutival.com)  
Email: [info@dr365pharmaceutival.com](mailto:info@dr365pharmaceutival.com)  
[dr365pharmaceutival@gmail.com](mailto:dr365pharmaceutival@gmail.com)  
Mob: 9594-70-2789 / 77-000-100-03

**BRAND :**  
**GE BrightSpeed 16 CT**

16 slice CT Scan System with Standard Accessories Year of Manufacture: December, 2012  
GE Bright Speed 16 with New Tube changed (Performix) in April, 2020 with Tube: 124,918 scan seconds as of 13th January, 2021



**HEAD OFFICE**



System Health Information						
Item	Information			Status		
Gantry Revolutions	4087403			Check PM Schedule		
Slipring Communication Bit Error Rate	0			OK		
DIP Errors	985e-15			OK		
Gantry Balance	0.25 kg-m			Balanced		
	0.10 kg-m2					
Item	Past Five Days			Today		
	Avg.	Min.	Max.			
Tube Spits / 10 Scan Seconds	0.5	0.05	1.5	Spins: 32, Scan Seconds: 593.3		
Slipring Brush Disconnects	0.00	0	10			
Item	Past Five Days			Today		
	Avg.	Min.	Max.	Avg.	Min.	Max.
Detector Temperatures	36.07	35.83	36.30	36.07	35.83	36.30
DAS ADB Temperatures	34.91	32.37	36.47	34.75	32.68	36.12
DAS IFB Temperatures	38.40	34.68	42.93	38.34	34.87	42.12
Gantry Ambient Temperatures	34.47	32.53	33.30	34.65	32.87	36.47

Auto Update  Update  
Frequency (sec) 10

Status: Successful with Display: 12'

Product Software Option Keys:

Smart Prep	Permanent
Patient-16-slice	Permanent
AutoA	Permanent
Direct-3D	Permanent
Direct-MPR	Permanent
Helical Tilt	Permanent
3000 Image Series	Permanent
Data Export	Permanent
CopyComposer	Permanent
VolumeViewer	Permanent
NeuroFilter	Permanent
CTColonPro	Permanent
CTPerfusion4Neuro	Permanent
AutoBone	Permanent
Connect Pro	Permanent
---End Of List---	